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WWW.KARMANIAPARSGEN.IR

AChE Activity Assay Kit
Catalog Nos. KPG-SOD

Acetylcholinesterase (AChE) is a cholinergic enzyme found mainly in postsynaptic neuromuscular junctions, especially in muscles and nerves. This enzyme rapidly breaks down or hydrolyzes acetylcholine (a natural neurotransmitter) into acetic acid and choline. This enzyme is also found in serum and red blood cells.

Contents of the kit:

1. Acetylcholinesterase reagent (CN. KPG-AChER)
2. Standard (CN. KPG-AchES)

Required items that are not included in the kit:

1. Sampler and sampler head, 2. Centrifuge, 3. ELISA reader with a wavelength of 560 nm, 4. Distilled water, 5. 96% ethanol.

Preparation of solutions:

To prepare the acetylcholinesterase reagent, mix 3 ml of 96% ethanol with 7 ml of distilled water and add it to the container of the acetylcholinesterase reagent and mix well. Be careful, if the reagent does not dissolve well (due to the difference in the quality of ethanol), try to dissolve the reagent using an ultrasonic bath or a sonicator probe. If the mentioned devices are not available, vortex well (at least 3 minutes).

The prepared reagent is stable for 72 hours in a dark place.

Preparation of samples

If the aim is to check the activity of total acetylcholinesterase, dissolve 10 microliters of whole blood containing anticoagulant with 8 ml of distilled water and after vigorous mixing, incubate for 5 minutes at room temperature.

If the aim is to investigate the activity of acetylcholinesterase in pure RBC, first wash 1 cc of blood with 10 cc of physiological serum for three times, then dissolve 100 microliters of the washed blood with 8 ml of distilled water and after Mix vigorously, incubate for 5 minutes at room temperature. Next, dilute the blood at a ratio of 1/20 with PBS buffer.

To check the activity of acetylcholinesterase in the tissue sample, first homogenize 50 mg of the tissue in 2 ml of PBS buffer, then centrifuge for 5 minutes at 5000 RPM and use the supernatant for testing.

Protocol:

1. Keep the kit components at room temperature for 20 minutes before starting the test.
2. Add 100 μ L of sample (prepared whole blood, washed or tissue supernatant) to a clean, sterile microtube.
3. Add 100, 50, 25, 12, 6 and 3 microliters of the standard to separate vials and name the standards 1 (72 U/mL), 2 (36 U/mL), 3 (U/mL 18), name 4 (9 U/mL), 5 (4.5 U/mL), and 6 (2.25 U/mL).
4. Add 100 microliters of distilled water to the blank vial.
5. Add 100 μ L of acetylcholinesterase reagent to all microtubes and incubate for 40 minutes at 37°C.
6. Add 100 microliters of colored produced solution to the ELISA plates and measure at a wavelength of 560 nm.

First, decrease the OD values of the blank from the samples and standards, then draw a graph and obtain the values of the samples based on the graph. Please note that the normal values of acetylcholinesterase activity are 18-8 U/mL.

Safety instructions: The solutions used in this kit are dangerous for human tissue, so be sure to work with gloves and gloves, and in case of contact with skin, eyes, etc., wash with plenty of water and go to the hospital for additional treatment.

Similar kits in the field of oxidants and antioxidants produced by Karmania Pars Gene:

1. Total Antioxidant Activity Test Kit (FRAP)
2. Catalase activity testing kit
3. MDA activity review kit
4. SOD activity test kit
5. Paraonase-1 activity testing kit
6. Protein carbonyl testing kit

Karmania Pars Gene produces kits:

1. ELISA to measure cytokines
2. Extraction of DNA and RNA by column, sediment and magnetic nanoparticle methods
3. Measurement of miRNAs
4. Products related to cell culture such as FBS, culture medium and cancer and mesenchymal stem cells
5. Microtubes, Real-Time PCR devices
6. DNA and RNA purification columns
7. Disposable laboratory products such as urine bottle and...

Contact us

00989132926113-00982191692296

karmaniaparsgene@gmail.com

WWW.KARMANIAPARSGEN.IR