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۰۹۱۳۲۹۲۶۱۱۳-۰۳۴۳۴۲۰۸۰۲۵

WWW.KARMANIAPARSGEN.IR

Ferric Reducing Ability of Plasma (FRAP) assay kit
Catalog Nos. KPG-FRAP

Contents of the kit:

1. Solution A (CN. KPG-Sa)
2. Solution B (CN. KPG-Sb)
3. Solution C (CN. KPG-Sc)
4. Solution D (CN. KPG-Sd)
5. Solution E (CN: KPG-Se)
6. Standard (CN. KPG-StdFRAP)

Required items that are not included in the kit:

1. Sampler and sampler head
2. Centrifuge
3. ELISA reader containing a wavelength of 593 nm

Preparation of solutions:

- Add 2.5 ml of solution E to the bottle of solution A and mix well.
- To prepare solution B, add 2.5 mL of deionized distilled water (DDW) to the bottle of solution B and mix well.
- To prepare solution C, add 15 ml of Solution D to the bottle of solution C and mix well.

The prepared solutions A, B and C are stable for 3 months in a dark place.

Preparation of working solution:

Solutions A, B and C should be mixed immediately before the test in a ratio of 1/10 to create a working solution. For example, if you need to examine 2 samples, prepare 300 μ L of working solution as follows: mix 25 μ L of solution A, 25 μ L of solution B and 250 μ L of solution C. The mixture is only stable for 30 minutes and should be prepared immediately before testing.

Preparation of standards:

Add 25 mL of deionized distilled water (DDW) to the standard bottle and mix well to obtain a 10 mM Fe²⁺ solution. Mix 2 mL of the standard containing 10 mM Fe²⁺ with 8 mL of distilled water to create a 2 mM Fe²⁺ standard, and then make serial dilutions using the table below to create standards 1 to 6 using a 2 mM solution. create:


Standards are stable at room temperature for 3 months in the dark.

Std	Concentrations (Fe ²⁺)	H ₂ O	Fe ²⁺ 2mm
1	1 mM (1000 μ M)	500 μ L	500 μ L
2	0.7 mM (700 μ M)	650 μ L	350 μ L
3	0.5 mM (500 μ M)	750 μ L	250 μ L
4	0.3 mM (300 μ M)	850 μ L	150 μ L
5	0.2 mM (200 μ M)	900 μ L	100 μ L
6	0.1 mM (100 μ M)	950 μ L	50 μ L

Test method:

1. Keep the kit components at room temperature for 20 minutes before starting the test.
2. Add 10 microliters of sample, standards and distilled water respectively to the wells of the ELISA plate corresponding to the sample, standards and blank.
3. Add 140 μ L of working solution to all wells and mix for 30 seconds.

Incubate the plate at 37°C for 5 minutes and then measure the OD of the samples along with the standards and the blank at a wavelength of 593 nm.

Sample standard curve. Please draw the standard curve yourself and do not use the curve.

Safety instructions:

The solutions used in this kit are hazardous to human tissue, therefore, in case of contact with skin, eyes, etc., wash with water and refer to hospital for additional treatment.

Contact us

00989132926113-00982191692296

karmaniaparsgene@gmail.com

WWW.KARMANIAPARSGEN.IR