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MDA Assay Kit (Catalog Nos. KPG-Nob)

Contents of the kit:

1. Solution A (CN. KPG-MDAa)
 2. Solution B (CN. KPG-MDAb)
 3. Solution C (CN. KPG-MDAc)
 4. Standard (CN. KPG-MDAs)
 5. Required items that are not included in the kit:
1. Sampler and sampler head, 2. Centrifuge, 3. ELISA reader with a wavelength of 535 nm

Preparation of solutions:

- Add 10 cc of distilled water to solution A and mix well.
- Add 10 cc of distilled water to solution B and mix well.

The composition of solution A and B is stable for one month after adding water.

Preparation of working solution:

Prepare the working solution as needed. To prepare the working solution, mix the A/B/C solution in a ratio of 2/1/1. For example, mix 2 cc of solution A with 1 cc of solution B and 1 cc of solution C.

The working solution has a stability of only one hour.

Preparation of standard solution:

To prepare the standard solution, dilute 20.5 microliters of the standard solution in the kit with 1 mL of 96% ethanol and then add to 49 mL of distilled water. The current solution contains 500 μ mol of standard. For the main standard, you must prepare a 50 micromol solution, that is, you must dilute the obtained standard solution to a ratio of 1/10. For this purpose, dissolve 200 microliters of the standard with 800 microliters of distilled water. Prepare subsequent standards at dilutions of 10, 5, 2.5, 1.25, 0.62, 0.31 and 0.15 μ mol.


Test method:

1. Add 200 microliters of the working solution to all microtubes related to standards, blanks, and samples.
2. Add 100 microliters of standards and 100 microliters of distilled water to the negative control and 100 microliters of the sample to the corresponding vial.
3. Incubate the microtubes for 60 minutes at 90°C.
4. Incubate for 10 minutes at room temperature until the temperature of the samples reaches room temperature.
5. Centrifuge the microtubes at 1000 RPM for 10 minutes.
6. Add 150 μ L of supernatant of standards, blanks, and samples to ELISA wells and immediately read at 535 nm.

Interpretation of results

After collecting the OD of the standards, plot the curve as shown in the standard sample and calculate the sample value using the standard curve and the formula $y=ax+b$.

To calculate the MDA concentration as μ M/mg pr or nmol/mg pr, please divide the MDA amount by the protein concentration. You can use Total Protein assay kit from Karmania Pars Gene Company.

An example of a standard diagram of the MDA activity monitoring kit.

Safety instructions: The solutions used in this kit are dangerous for human tissue, so be sure to work with gloves and gloves, and in case of contact with skin, eyes, etc., wash with plenty of water and go to the hospital for additional treatment.

Similar kits in the field of oxidants and antioxidants produced by Karmania Pars Gene:

1. Total Antioxidant Activity Test Kit (FRAP)
2. Catalase activity testing kit
3. Kit to check the amount of NO
4. SOD activity test kit
5. Paraonase-1 activity testing kit
6. Protein carbonyl testing kit

Karmania Pars Gene produces kits:

1. ELISA to measure cytokines
2. Extraction of DNA and RNA by column, sediment and magnetic nanoparticle methods
3. Measurement of miRNAs
4. Products related to cell culture such as FBS, culture medium and cancer and mesenchymal stem cells
5. Microtubes, Real-Time PCR devices
6. DNA and RNA purification columns
7. Disposable laboratory products such as urine bottle and...

Contact us

00989132926113-00982191692296

karmaniaparsgene@gmail.com

WWW.KARMANIAPARSGEN.IR