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Required items that are not included in the kit:

1. **Sampler and sampler head,**
2. **Centrifuge,**
3. **ELISA reader containing 450 nm wavelength**

Preparation of solutions:

Add solution A completely to solution B.

To prepare solution D, add one mL of distilled water to the vial and mix well.

The combination of solution A and B is stable for up to 72 hours and solution D is stable for one hour after adding distilled water.

Preparation of standard solution:

To prepare the standard solution, mix D and C in a ratio of 1/1. For example, mix 100 µl of solution D and 100 µl of solution C for each test.

Protocol:

1. Keep the kit components at room temperature for 20 minutes before starting the test.
2. To make a standard with dilutions of 100, 50, 25, 12.5 and 6.25 micromol, add 30, 15, 7.5, 3.75 and 1.85 microliters of the prepared mixture to standard microtubes 1 to 5, respectively.
3. Add 100 microliters of the mixture of solutions A and B to all the microtubes, including the standard, blank, and samples (according to the table below).
4. Add 30 microliters of sample to the microtube corresponding to the samples and 30 microliters of distilled water (DW) to the blank (according to the table below).

	S1	S2	S3	S4	S5	Sample	DW
A and B mix	100 µL	100 µL	100 µL	100 µL	100 µL	100 µL	100 µL
Prepared standard	30 µL	15 µL	7.5 µL	3.75 µL	1.85 µL	-	-
Test	-	-	-	-	-	30 µL	-
Blank	-	-	-	-	-	-	30 µL

5. Mix the microtubes vigorously and incubate for 10 minutes at room temperature in the dark.
6. Add 150 microliters of solution E to all microtubes and mix.
7. Add 150 µL of standards, blanks, and samples to ELISA wells and immediately read at 450 nm.

Safety instructions: The solutions used in this kit are dangerous for human tissue, so be sure to work with gloves and gloves, and in case of contact with skin, eyes, etc., wash with plenty of water and go to the hospital for additional treatment.

Similar kits in the field of oxidants and antioxidants produced by Karmania Pars Gene:

1. Total Antioxidant Activity Test Kit (FRAP)
2. Catalase activity testing kit
3. MDA activity review kit
4. SOD activity test kit
5. Paraonase-1 activity testing kit
6. Protein carbonyl testing kit

Karmania Pars Gene produces kits:

1. ELISA to measure cytokines
2. Extraction of DNA and RNA by column, sediment and magnetic nanoparticle methods
3. Measurement of miRNAs
4. Products related to cell culture such as FBS, culture medium and cancer and mesenchymal stem cells
5. Microtubes, Real-Time PCR devices
6. DNA and RNA purification columns
7. Disposable laboratory products such as urine bottle and...

Contact us

00989132926113-00982191692296

karmaniaparsgene@gmail.com

WWW.KARMANIAPARSGEN.IR