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WWW.KARMANIAPARSGEN.IR

Total protein Assay Kit Catalog Nos. KPG-ToP

The BCA protein quantification kit is a simple method to determine the concentration of total proteins. In this kit, Cu^{2+} changes color by binding to proteins and creates an intensely colored complex with bicinchoninic acid, which is read at a wavelength of 562 nm.

Contents of the kit:

1. Solution A (CN. KPG-TOPA)
2. Solution B (CN. KPG-ToPB)
3. Standard (CN. KPG-ToPS)

Required items that are not included in the kit:

1. Sampler and sampler head
2. Centrifuge
3. Eliza plate
4. ELISA reader containing 562 nm wavelength

Preparation of working solution:

- To prepare the working solution, mix solution A with solution B at a ratio of 1.50. For example, mix 10 microliters of solution A with 490 microliters of solution B.
- The prepared working solution is stable for 24 hours in a dark place.

Preparation of standards:

First, prepare six sterile microtubes and label them with standard names 1, 2, 3, 4, 5 and 6. Add 450 microliters to the first microtube and 250 microliters to the remaining microtubes of the buffer in which your protein is dissolved.

Next, add 50 microliters of protein to the first microtube and mix well to prepare a 1000 $\mu\text{g/mL}$ standard. Next, add 250 microliters of the first standard (1000 $\mu\text{g/mL}$) to the second standard and mix well to prepare the second standard with a concentration of 500 $\mu\text{g/mL}$. Do the same until the final standard, until the concentrations of μg 250, 125, 62.5 and 31.25/mL should be prepared.

Protocol:

1. Keep the kit components at room temperature for 20 minutes before starting the test.
 2. Add 25 μL of sample and standards to separate microtubes and 25 μL of distilled water as a control to separate microtube.
 3. Add 225 microliters of working solution to all microtubes and incubate at 60°C for 60 minutes.
 4. Add 100 μL of the solution to all wells and measure at 562 nm wavelength. Pay attention that if you don't have 562 wavelength, you can also read at 545 nm wavelength.
- To calculate the amount of protein in the sample, first subtract the OD of the blank from the OD of the samples and standards, and then calculate the amount of protein by drawing the standard graph.

Safety instructions: The solutions used in this kit are dangerous for human tissue, so be sure to work with gloves and gloves, and in case of contact with skin, eyes, etc., wash with plenty of water and go to the hospital for additional treatment.

Similar kits in the field of oxidants and antioxidants produced by Karmania Pars Gene:

1. Total Antioxidant Activity Test Kit (FRAP)
2. Catalase activity testing kit
3. Kit to check the amount of NO
4. FRAP Assay test kit
5. Paraoxonase-1 activity testing kit
6. Protein carbonyl testing kit

Karmania Pars Gene produces kits:

1. ELISA to measure cytokines
2. Extraction of DNA and RNA by column, sediment and magnetic nanoparticle methods
3. Measurement of miRNAs
4. Products related to cell culture such as FBS, culture medium and cancer and mesenchymal stem cells
5. Microtubes, Real-Time PCR devices
6. DNA and RNA purification columns
7. Disposable laboratory products such as urine bottle and...

Contact us

00989132926113-00982191692296

karmaniaparsgene@gmail.com

WWW.KARMANIAPARSGEN.IR