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Plasmid extraction kit (CN: KPG-PLS)

A brief explanation of the Plasmid extraction kit:

To extraction of plasmid it is needed to exclude bacterial DNA and RNA. To the elimination, Solution II is used. Please perform the incubation accurately.

Kit component

Item	Catalogue number	volume
Solution I	KPG-SI	10 mL
Solution II	KPG-SII	15 mL
Solution III	KPG-SIII	20 mL
Solution IV	KPG-SIV	10 mL
Precipitation solution	KPG-PS	25 mL
Washing buffer I	KPG-WB	25 mL
DNase/RNase free water	KPG-DW	5 mL
High absorbance column tube	KPG-Column	50 pcs
Collection Tube	KPG-CT	100 pcs

Material needs

Distilled Water	Sampler tips
Centrifuge	Samplers

KPG-Plasmid kit was designed to extract and purify Plasmid from Bacteria (Gram-positive and negative bacteria). **KPG-Plasmid** kit is fast and easy to process and has a High column binding capacity. All reagents are stable at room temperature for one year. For extracting and purifying Plasmid from bacteria, culture bacteria to reach to Log phase, preferably. Extract and purify the plasmid as the following procedure.

Experimental procedure

1. Prepare a high-density bacterial suspension (at least 3×10^9 CFU/mL) from fresh colonies or broth culture in distilled water into a 2 mL Micro-tube. Vortex bacterial suspension for 5 s and centrifuge at $8000 \times g$ / 2 min
2. Throughout supernatant and suspend pellet in 200 μ L **S I solution**. Vortex at high speed and incubate at room temperature for 5 min.
3. Add 250 μ L **S II solution** and mix gently. Incubate at room temperature for a maximum of 1 min: **Attention; do not exceed the time over 1 min. It will be damaged to plasmid DNA.**
4. Add 350 μ L of **S III solution** to the tube and mix Micro-tube context by inverting (30 s).
5. Add 200 μ L **S IV buffer** and mix by inverting and then Centrifuge the Micro-tube at $13000 \times g$ / 5 min at room temperature.
6. Transfer carefully the supernatant to the new 2 mL micro-tube. Do not transfer the pellet.
7. Add 500 of **precipitation buffer** to the 2 mL micro-tube and mix by inverting.
8. Transfer the lysate to the 2 mL collection tube with a spin column tube and Centrifuge at $12000 \times g$ / 1 min at room temperature. Discard the 2 mL micro-tube and transfer the spin column tube into the new 2 mL collection tube. **"If the lysate remains, transfer it to the column and centrifuge again"**
9. Add 500 μ L **Washing buffer** and centrifuge at $12000 \times g$ / 1 min at room temperature.
10. Transfer the spin column tube to the DNase-free 1.5 mL Micro-tube
11. Add 50 μ L DNase-free water and incubate in a water bath at 65°C / 4
12. Centrifuge the Microtube at $12000 \times g$ / 2 min at room temperature
Discard the column tube, the pellet contains the DNA of the plasmid. Keep plasmid DNA at -20°C for later use.

Safety while using the kit

The solutions used in the kit have oxidizing and acidic properties. Avoid direct contact with skin and eyes. In case of contact with the mentioned tissues, wash with plenty of water and go to the nearest medical center.

Factory Address

Refsanjan, 20 km from Rafsanjan road to Kerman, Rafsanjan special economic zone, food and drug area
09132926113
09135025983
Karmaniaparsgene@gmail.com
www.karmaniaparsgen.ir

Other DNA extraction kits and RNA produced by Karmania Pars Gen:

All the following DNA extraction kits are available in all three column, sedimentary and magnet nanoparticles:

1. DNA extraction kits from blood and tissue by column
2. DNA extraction kits from gram-positive and negative bacteria
3. DNA extraction kit from mycobacterium
4. DNA extraction kit from Fungi
5. DNA extraction kit from virus
6. DNA extraction kit from HPV virus
7. DNA extraction kits from plant tissue

RNA extraction kits:

All the following RNA extraction kits are available in both column and sedimentary methods:

1. RNA extraction kits from blood and tissue
2. RNA extraction kits from gram-positive and negative bacteria
3. RNA extraction kit from mycobacterium
4. RNA extraction kit from Fungi
5. The RNA extraction kit from the virus
6. RNA extraction kits from plant tissue