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DNA extraction kit from stool

(CN: KPG-DNA-Stool)

A brief explanation of the DNA extraction kit from the stool:

To discharge DNA from biological samples, samples of waste must be discharged, followed by a completely homogeneous and subsequent cell lysis. There are many waste in the stool that is separated from the process in the first thigh. Using the enzymes of Lysozyme and the required buffers of all cells, including gram -negative and positive bacteria, it is completely lightened and DNA will attach to these filters by passing through the columns containing Silica and after treating with DNase/RNase Free water, DNA sample is separated from the filter.

Kit component

Item	Catalogue number	volume
Lysozyme	KPG-L	1 mL
KPG PreLysis Solution	KPG- KPS	10 mL
KPG Lysis solution	KPG-KLS	25 mL
Precipitation solution	KPG-PS	15 mL
Carrier	KPG-CR	1mL
Washing buffer	KPG-WB	25 mL
DNase free water	KPG-DW	5 mL
Filtered column tube	KPG-FColumn	50 pcs
High absorbance column tube	KPG-Column	50 pcs
Collection Tube	KPG-CT	100 pcs

Material needs

Distilled Water	Sampler tips
Centrifuge	Samplers

Sample Preparation Steps: To extract DNA, mix about 0.25 grams of homogenized stools in 800 microliters of distilled water and mix with a high-speed for 30 seconds, then make a quick rotation with microfuge (1-2 seconds) to deposit large particles. Transfer the fecal fluid to the filter column tube (KPG-FColumn) and centrifuge at 7000 rpm for 2 minutes. Collect the fluid passing through the filter to a 2ml sterile microtube and add 100 microliters of prelysis buffer (KPG-KPS). Next, add 20 microliters of Lysozyme (40 mg / ml) and mix for 10 seconds and then incubate at least 60 minutes at 37 ° C.

How to work with the kit to extract DNA:

1. Add 400 microliter KPG Lysis Buffer (KPG-KLS) to the previous step microtube and vortex it for 10 seconds.
2. incubate the suspension at 65 ° C for 5 minutes until the bacterial cells are fully lysed.
3. Centrifuge the microtube at a speed of 7,000 rpm 3 minutes and transfer the supernatant to 1.5 ml of sterile microtube.
4. Add 10 µL of Carrier (KPG-CR) and vortex for 5 seconds.
5. Add 300 microliter Precipitation Solution (KPG-PP) and vortex for 5 seconds.
6. Transfer the entire the tube to the DNA extraction column (KPG-Column) and centrifuge at 8000 rpm for 1 minute.
7. Transfer the column to the new collection tube and add 500 µm from the washing buffer (KPG-WB) to the column and centrifuge at 8000 rpm for 1 minute.
8. Discharge the contents of the collection tube and re -insert the column inside the collection tube and centrifuge at 12,000 rpm for 1 minute without adding any solution.
9. Transfer the column to the new collection tube and add 50 µL of distilled water and incubate for 4 minutes at room temperature. It is advisable to pre-warm distilled water to 65 ° C.
10. Centrifuge the column at 12000 for 1 minute. The solution passed through the column contains DNA.

Note: To check the accuracy and amount of DNA drainage, be sure to check the DNA OD at 260 and 280 nm wavelengths, and ensure that it is accurate by electrophoresis.

Safety while using the kit

The solutions used in the kit have oxidizing and acidic properties. Avoid direct contact with skin and eyes. In case of contact with the mentioned tissues, wash with plenty of water and go to the nearest medical center.

Factory Address

Refsanjan, 20 km from Rafsanjan road to Kerman, Rafsanjan special economic zone, food and drug area
09132926113
09135025983
Karmaniaparsgene@gmail.com
www.karmaniaparsgen.ir

Other DNA extraction kits and RNA produced by Karmania Pars Gen:

All the following DNA extraction kits are available in all three column, sedimentary and magnet nanoparticles:

1. DNA extraction kits from blood and tissue by column
2. DNA extraction kits from gram-positive and negative bacteria
3. DNA extraction kit from mycobacterium
4. DNA extraction kit from Fungi
5. DNA extraction kit from virus
6. DNA extraction kit from HPV virus
7. DNA extraction kits from plant tissue

RNA extraction kits:

All the following RNA extraction kits are available in both column and sedimentary methods:

1. RNA extraction kits from blood and tissue
2. RNA extraction kits from gram-positive and negative bacteria
3. RNA extraction kit from mycobacterium
4. RNA extraction kit from Fungi
5. The RNA extraction kit from the virus
6. RNA extraction kits from plant tissue

One of the most efficiently discharged kits: **Plasmid extraction Kit**